



REGULAR (adult) sized mask



	4	130mm x 53mm
then this... etc	5	130mm x 40mm
mount this first	6	130mm x 11mm
	7	130mm x 13mm
	8	130mm x 13mm
	9	130mm x 13mm
	10	130mm x 70mm



Print on card as thick as you think your printer can handle.

Cut out the parts 1-10 (previous pages) and cut AROUND the Sewing Jig (above).

Mount parts 1-10 onto the Sewing Jig with clear Sellotape so that the strips are hinged on the axis indicated. Use tape on both sides of each strip so they do not dislocate.

... Watch the video for "how to"

<https://www.daydreamfoundry.com/maskproduction/video>

Please read through the info found here before you start mass producing:

<https://www.daydreamfoundry.com/maskproduction/>

The choice of metal wire gauge is particularly important to get a good lasting fit around the nose. Stiffer is better than sloppy.

**Length of metal strip/wire 100mm**

# Material cutting template

REGULAR (adult) sized mask

175mm x 200mm

NB If using a vacuum bags as your filter material - Place the fabric in the Jig with the Bag inside facing down. Ie. the Bag Outside becomes the Mask Inside.